HARTMANN

# HydroTerapie

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
[1] Atkin, L. and Ousey, K. (2016). Wound bed preparation: A novel approach using HydroTherapy. *British Journal of Community Nursing* 21 (Suppl. 12), pp. S23-S28. [2] Ousey, K. et al. (2016). HydroTherapy Made Easy. *Wounds UK* 12(4). [3] Humbert, P. et al. (2014). Protease-modulating polyacrylate-based hydrogel stimulates wound bed preparation in venous leg ulcers – a randomized controlled trial. *Journal of the European Academy of Dermatology and Venereology* 28(12), pp. 1742-1750. [4] Smola, H. (2016). Simplified treatment options require high-performance dressings – from molecular mechanisms to intelligent dressing choices. EWMA 2016. Bremen, 11-13 May, 2016. [5] Smola, H. et al. (2016). Hydrated polyurethane polymers to increase growth factor bioavailability in wound healing. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing*. London, 3 March, 2016. [6] Eming, S., Smola, H., Hartmann, B. et al. (2008). The inhibition of matrix metalloproteinase activity in chronic wounds by a polyacrylate superabsorbent. *Biomaterials* 29: 2, 2932-2940 [7] ISBN Monograph. 978-1-944788-63-6-ISBN Services [8] Kaspar, D. (2011). Therapeutic effectiveness, compatibility and handling in the daily routine of hospitals or physicians's practices. HARTMANN Data on file: Hydro-Responsive Wound Dressing (HRWD) and AquaClear Technology are trademarks of HARTMANN [9] Ousey, K. et al. (2016). Hydro-Responsive Wound Dressings simplify T.I.M.E. wound management framework. *British Journal of Community Nursing* 21 (Suppl. 12), pp. S39-S49. [10] Spruce, P. et al. (2016). Introducing HydroClean® plus for wound-bed preparation: a case series. *Wounds International* 7(1), pp. 26-32. [11] Ousey, K. et al. (2016). HydroClean® plus: a new perspective to wound cleansing and debridement. *Wounds UK* 12(1), pp. 94-104. [12] Atkin, L. and Rippon, M. (2016). Autolysis: mechanisms of action in the removal of devitalised tissue. *British Journal of Nursing* 25(20), pp. S40-S47. [13] Kaspar, D. et al. (2015). Economic benefit of a polyacrylate-based hydrogel compared to an amorphous hydrogel in wound bed preparation of venous leg ulcers. *Chronic Wound Care Management and Research* 2, pp. 63-70. [14] Bullough, L. et al. (2016). A multi-centre 15 patient evaluation of a Hydro-Responsive Wound Dressing (HRWD) - HydroClean® plus. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing*. London, 3 March, 2016. [15] Scherer, R. et al. (2015). HydroTherapy®. Application study. Heidenheim: Paul Hartmann AG. [Data on file]. [16] Chadwick, P. and Haycocks, S. (2016). The use of Hydro-Responsive Wound Dressing for wound bed preparation in patients with diabetes. *Wounds UK Annual Conference*. Harrogate, 14-16 November, 2016. [17] Knowles, D. et al. (2016). HydroTherapy® wound healing of a post amputation site. *Wounds UK Annual Conference*. Harrogate, 14-16 November, 2016. [18] Zollinger, C. et al. (2014). HydroTherapy®. Application Study. Heidenheim: Paul Hartmann AG. [Data on file]. [19] Colegrave, M. et al. (2016). The effect of Ringer's solution within a dressing to elicit pain relief. *Journal of Wound Care* 25(4), pp. 184-190. [20] O'Brien, D. and Clarke, Z. (2016). The patient experience with a Hydro-Responsive Wound Dressing (HRWD) – HydroClean® plus. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing*. London, 3 March, 2016. [21] Jones, T. and McCracken, K. (2016). HydroClean® plus assists healing of leg ulcers for a patient with systemic lupus erythematosus. *Wounds UK Annual Conference*. Harrogate, 14-16 November, 2016. [22] Kaspar, D. et al. (2008). Efficacité clinique du pansement irrigo-absorbant HydroClean® active contenant du polyacrylate superabsorbant dans le traitement des plaies chroniques – étude observationnelle conduite sur 221 patients. *Journal des Plaies et Cicatrisations* 13(63), pp. 21-24. [24] Smola, H. et al. (2014). Hydrated polyurethane polymers to increase growth factor bioavailability in wound healing. *EORS Congress*. Nantes, 2-4 July, 2014. [25] Spruce, P. and Bullough, L. (2016). HydroTac®: case studies of use. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing*. London, 3 March, 2016. [26] Spruce, P. et al. (2016). A case study series evaluation of HydroTac®. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing*. London, 3 March, 2016. [27] Smola, H. et al. (2016). From material science to clinical application – a novel foam dressing for the treatment of granulating wounds. *HydroTherapy Symposium: A New Perspective on Wound Cleansing, Debridement and Healing* [28] Research Results Product Test. Data on file. [29] Pharmacy prescription data 2015. [30] www.vidal.fr/parapharmacie/irriter/id9233/hydrotac\_transparent\_pans\_9eL\_non\_adhesif/.21-02-14. [31] Data on file. [32] Knestele, M (2004) The treatment of problematic wounds with HydroClean advance - tried and tested over many years in clinical practice. Údaje HARTMANN ve spisu. [33] Zöllner, P (2010). Observational study Ideal wound environment during granulation and epithelization achieved with HydroTac. Údaje HARTMANN ve spisu. [34] Meuleneire F (2011b) Clinical evaluation of TenderWet plus. Dokument HARTMANN [Údaje ve spisu]. [35] Bruggisser, R. (2005). Bacterial and fungal absorption properties of a hydrogel dressing with a superabsorbent polymer core. *J Wound Care* 14, 438-42. [36] Ellermann, J (2015). HydroClean 2.0: Design validation customer/user interviews. Internal Report, International Marketing Department. [37] Data on file: 27. Z+SilBor-der\_benchmark. [38] Call, E., et al., Comparing fluid handling and microclimate conditions under Superabsorbent polymer and superabsorbent foam dressings over an artificial wound. *WCET Journal*, 2019.

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